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



TAP <ファイル>


TAP <オプション>

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Patient information and driving unit settings

Patient ID	EXCOR Active	Language for driving unit	English
Patient weight	36.5 kg	Child lock	Off
Implantation date	07 May 2024	Time of self-test	10:22 AM

Patient information and driving unit settings > Blood pumps and cannulae

● ○ ○

Logged in as expert

L 

 R 

EXCOR Active

10:54 AM
07 May 2024

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Patient weight

2.0 kg

255.0 kg

20.6 kg

Apply

Cancel

Logged in as expert

  L  R 

EXCOR Active

10:53 AM
07 May 2024

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Patient weight

2.0 kg 255.0 kg

20.6 kg


1 2 3
4 5 6
7 8 9
C 0 .

Apply Cancel

Logged in as expert L R EXCOR Active 10:53 AM 07 May 2024

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The screenshot displays the 'Patient information and driving unit settings' screen. At the top, a navigation bar contains icons for menu, data transfer, patient status, pump settings, and alerts. The main content area is titled 'Patient information and driving unit settings' and contains six input fields arranged in two columns. The left column includes 'Patient ID' (BH1-087), 'Patient weight' (20.0 kg), and 'Implantation date' (05 Feb 2025). The right column includes 'Language for driving unit' (English), 'Child lock' (Off), and 'Time of self-test' (09:19 AM). A hand icon with a blue arrow points to the 'Time of self-test' field. Below the input fields, a progress indicator shows three dots, with the first dot filled, indicating the current step. To the right of the progress indicator is a link '> Blood pumps and cannulae'. At the bottom, a status bar shows 'Logged in as expert', pump status icons for L and R, the patient ID 'BH1-087', and the time '09:23 AM' on '05 Feb 2025'.

Field	Value
Patient ID	BH1-087
Patient weight	20.0 kg
Implantation date	05 Feb 2025
Language for driving unit	English
Child lock	Off
Time of self-test	09:19 AM


Patient information and driving unit settings > Blood pumps and cannulae

Logged in as expert

BH1-087 09:23 AM 05 Feb 2025

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The screenshot displays the Berlin Heart panel PC interface. At the top, a navigation bar contains several icons: a menu icon, a pump icon, a cannula icon, a settings icon (highlighted with a white triangle), and a bell icon. The main title 'Blood pumps and cannulae' is centered below the navigation bar. The interface shows three main sections: 'Operating mode' set to 'BVAD', 'Blood pump [mL]' with values '30' and '25', and 'Cannula inflow / outflow [mm]' with values '9 | 9'. At the bottom, a status bar includes 'Patient information and driving unit settings <', 'Blood pumps and cannulae' (with a selected dot), '> Settings' (with a hand cursor pointing to it), and a chat icon. The bottom-most bar shows 'Logged in as expert', pump status icons, 'BH1-087', and the date/time '06:22 PM 06. Jan 2021'.

Operating mode

BVAD

Blood pump [mL]

30 | 25

Cannula inflow / outflow [mm]


9 | 9

Patient information and driving unit settings < Blood pumps and cannulae > Settings

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